

Ergonomic

The innovative system with laces, allows to tie the apron directly to the patient's neck.



Colourful

They come in numerous colours that can be matched with other Monoart® products.



Protective

Practical shape that perfectly covers the patient's neck and shoulders.

Monoart® Apron PG30



A roll of disposable aprons for patient protection. They are made of a layer of highly absorbent paper and a layer of waterproof polyethylene, which are bonded and embossed. Unlike traditional towels, they provide perfect protection for patients thanks to their large size and the specially shaped neck and shoulder sections.



Dimensions:

Roll of 80 aprons:
610x530 mm

Roll of 60 aprons:
810x530 mm



Colours:


Lime, green, emerald green, blue, light blue, lilac, pink, burgundy, orange, yellow, white, grey



Packaging:

Carton of 6 rolls

Instructions for use



Hygienic

No glue is used. Embossing and heating are used to bond the polyethylene film and fabric together.



Absorbent

It guarantees high speed of absorbency to protect the patient.