



Eco-friendly

Monoart® Aspirator Tips can go through up to 100 cycles in an autoclave at up to 134°C.



Comfortable

With a smaller contact surface thanks to the smaller end “brims” and fewer curves, the tips offer greater patient comfort.



Effective

The smooth interiors of the tips allow liquids to flow more easily.

Monoart[®] EM19 Evo Aspirator Tips



Monoart[®] EM19 Evo universal tips for surgical aspirators are recommended for procedures on adults. Ergonomic, non-slip grips make work easier for staff, while the linear shape keeps them quiet and promotes better hygiene and cleaning. EM19 Evo tips give outstanding patient comfort thanks to their rounded shape. They can be sterilized in autoclaves at up to 134°C.

**Classification:**

Medical Device Class IIa

**Dimensions:**

Connection diameter:
16 mm; Length: 114 mm

**Colours:**

Lime, green, blue, lilac,
pink, burgundy, orange,
yellow, grey

**Packaging:**

Box of 10 pieces
of the same colour

**Safe**

The grooves are more extensive and reach all of the way to the end, for better grip.

**Ergonomic**

The tips are shorter and bent towards the end, for better control of operations on patients.

**Silent**

The angle of approximately 16° in the tips makes aspiration quieter.

Monoart® EM21 Evo Aspirator Tips



Monoart® EM21 Evo universal tips for surgical aspirators are recommended for procedures on cavities that are difficult to access. Ergonomic, non-slip grips make work easier for staff, while the linear shape keeps them quiet and promotes better hygiene and cleaning. They can be sterilized in autoclaves at up to 134°C.



Classification:

Medical Device Class IIa



Dimensions:

Connection diameter:
16 mm; Length: 129 mm



Colours:

Lime, green, blue, lilac,
pink, burgundy, orange,
yellow, grey



Packaging:

Box of 10 pieces
of the same colour



Multisize

EM 19 Evo


Connection diameter: 16 mm
Length: 114 mm

EM 21 Evo

Connection diameter: 16 mm
Length: 129 mm

EM 40

Connection diameter: 11 mm
Length: 200 mm



Monoart[®] EM40 Aspirator Tip



Universal moulded tip for surgical aspirator.
Can be sterilised in autoclave at 134°C.



Classification:

Medical Device Class IIa



Dimensions:

Connection diameter:
11 mm; Length: 200 mm



Colour:

Grey



Packaging:

Box of 10 pieces

Monoart[®] Adapter for Aspirator Tip



Safe and practical adapter for EM40 aspirator tip.
Can be sterilised in autoclave at 134°C.



Classification:

Medical Device Class IIa



Packaging:

Bag of 10 pieces