



Safe

They provide an exceptional protective barrier against contamination and infection and they do not contain lubricant powders, chemical residue and endotoxins, which are often responsible for the respiratory symptoms caused by latex.



Certified

They comply with strict quality control criteria:

- AQL 1.5
- Level G1
- Low levels of protein residue (< 50 mg/g).

They are certified in accordance with the international directive 93/42/EEC, as amended by directive 2007/47/EC.



Durable

They are water, acid and alkali resistant and they provide good protection against micro-organisms and the majority of caustic substances and cleaning products.

Monoart® Disposable Latex Gloves



Safe, comfortable and durable: Monoart® disposable latex gloves provide hygienic, comfortable protection for hands. They are easy to wear, stretchy and guarantee outstanding tactile sensitivity. In addition, the rough micro-texture on the surface gives superb grip in all conditions, so they are ideal for working with complete peace of mind during operations of all kinds. Monoart® latex gloves come in a range of bright colours to give a unique look to every dental practice.

**Classification:**

Medical Devices Class I, pursuant to EN 455-1, EN 455-2, EN 455-3 and EN 455-4

**Weight:**

6.2 ± 0.3 gr

**Available sizes:**

XS, S, M, L

**Colours:**

Lime, blue, lilac, pink

**Packaging:**

Box of 100 gloves
Carton of 10 boxes

**Comfortable**

They offer outstanding tactile sensitivity and they are easy to put on thanks to the polymer lining and the folding cuffs.

**Free**

They come in bright colours to give a unique look to every dental practice. Enjoy peace of mind at work thanks to these safe, resistant and comfortable gloves.