



## Practical

The copper-clad wire is integrated in the tube to keep it at the set angle and ensure it won't protrude when bent.



## Safe

Made of non-toxic PVC without DOP and phthalates, Class IIa Medical Devices with no expiry date.



## Comfortable

With rounded fixed or removable tips that won't irritate mucous membranes.

# Monoart<sup>®</sup> Saliva Ejectors EM15



Disposable, flexible saliva ejectors made of non-toxic, transparent PVC. The integrated metal wire perfectly holds them at the required angle and never protrudes during the bending process. Fixed or removable rounded tips are ideal for work on irritated mucous membranes.

**Classification:**

Medical Device Class IIa

**Dimensions:**

Diameter: 6.5 mm  
Length: 12.5 or 15 cm

**Colours:**

Blue, green, lime, yellow, transparent, orange, pink, burgundy, lilac

**Packaging:**

Bag of 100 pieces  
Carton of 50 bags

# Monoart<sup>®</sup> Adapter for Saliva Ejector



Safe, practical and autoclavable adapter for saliva ejector. Can be sterilized in autoclave at 134°C.

**Classification:**

Medical Device Class IIa

**Packaging:**

Bag of 10 pieces